<b>(D)</b> Gebrauchsinformation .....	3
<b>(GB)</b> Instructions for Use .....	11
<b>(F)</b> Notice d'utilisation .....	18
<b>(E)</b> Instrucciones de uso .....	26
<b>(I)</b> Istruzioni per l'uso .....	34
<b>(P)</b> Instruções de uso .....	41
<b>(NL)</b> Gebruiksaanwijzing .....	48
<b>(S)</b> Bruksanvisning .....	56
<b>(DK)</b> Brugsanvisning .....	63
<b>(N)</b> Bruksveiledning .....	70
<b>(FIN)</b> Käyttöohje.....	77
<b>(GR)</b> Πληροφορίες χρήσης .....	84
<b>(PL)</b> Sposób użycia .....	92
<b>(CZ)</b> Návod k použití .....	100
<b>(HR)</b> Használati utasítás .....	107
<b>(LV)</b> Informācija lietotājam .....	115
<b>(LT)</b> Vartotojo informacija .....	122

### Produktbeschreibung

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Zubehör

**Description**

iBOND Self Etch is a light-curing self-priming one-component bonding agent for use in combination with adhesive restorations.

Separate conditioning (etching) of the enamel and dentine is not required (however, use of an additional etching gel on the enamel before application of iBOND Self Etch will not have a negative influence on the bond strength).

iBOND Self Etch was developed for bonding composite resin materials (e.g. composite, compomer, Polyglas<sup>®</sup>) to the hard tooth structure. iBOND Self Etch etches, primes, bonds and desensitises in one step.

For dental use only. Do not use for indications or applications that are not specifically noted in the instructions for use.

Read the instructions for use carefully before use!  
Follow the safety instructions when using iBOND Self Etch!

## Composition

iBOND Self Etch is an acetone/water-based formulation of light-activated methylacrylate resins.

## Indications

1. Bonding of direct light-cured composite restorations (including Polyglas® and compomers).
2. Bonding of indirect restorations in combination with a light-curing luting cement: porcelain, Polyglas® and composite restorations (inlays, onlays, veneers, crowns).
3. Sealing hypersensitive areas of teeth.

## Storage

Do not store above room temperature (23 °C) during daily use. If iBOND Self Etch is not yet opened (bottle version), it must be stored in a refrigerator at 4–10 °C (40–50 °F). Bring to room temperature and shake bottle briefly prior to use. Replace and tighten bottle cap, immediately after use. If not stored correctly the product may lose its effectiveness early. If one of the following is noted, it is a sign that the material is losing its bonding efficacy and should no longer be used:

- Nothing comes out of the bottle when squeezed although there is material left inside.
- Transparent solid or gel-like coagulation is detected in the liquid when dispensed.
- Threads are visible in the material during application with the micro brush.

## Application

### 1. Bonding of direct light-cured composite restorations (including Polyglas® and compomers).

#### 1.1. Preparation

Do not use iBOND Self Etch directly on exposed pulp tissue.

- Clean tooth with an oil and fluoride-free polishing paste.
- Prepare cavity for bonding.
- Rinse cavity with water and air dry.
- Isolate the cavity (the use of a rubber dam is strongly recommended).
- If indicated, place sub-base (e.g. glass ionomer cement).

#### 1.2. Dispensing

Bottle application:

- Shake bottle briefly before opening.
- Dispense iBOND Self Etch into the well (1 drop for small cavities, 2 drops for larger cavities).
- Replace bottle cap immediately after dispensing.
- Use iBOND Self Etch immediately after dispensing (no later than 3 minutes).

Single dose application:

- Make sure that the single dose is opened immediately prior to application.
- When removing the liquid rotate the applicator tip or brush briefly in the single dose.

### 1.3. Application

Roughen or grind enamel before application of iBOND Self Etch. Etching/bonding is less effective on unground enamel. On sclerotic dentine separate etching with phosphoric acid for 30 s is recommended.

- Immediately after dispensing apply a copious amount of iBOND Self Etch to the entire cavity surface and margin with the applicator tip or brush. Make sure that the cavity margin is also sufficiently coated with iBOND Self Etch liquid.
- After application agitate the adhesive slightly for 20 s. Agitation during the dwell time improves demineralisation and diffusion. Protect the bonding layer from contamination (e.g. by blood or saliva).
- Carefully air-dry iBOND Self Etch with a flow of oil-free air (may take 5–10 s or longer depending on the geometry of the cavity). The objective is to evaporate the solvent and water from the bonding layer without removing the active ingredients from the tooth surface.

A strong air stream at the start of air-drying will dilute the bonding agent and may result in poor adhesion.

- The surface must be visibly glossy, both after application of iBOND Self Etch and after evaporation of the solvent. The complete cavity surface must be fully coated. If the cavity surface does not appear shiny, apply iBOND Self Etch as described above a second time.
- Light-activate iBOND Self Etch for 20 s with a halogen or LED light-curing unit. The light-curing time assumes the use of a Heraeus Kulzer Translux® light-curing unit or other dental curing light of comparable intensity (min. 400–500 mW/cm<sup>2</sup>). If plasma light-curing units (with an output greater than 1200 mW/cm<sup>2</sup>) are used, the curing time can be reduced by 8 s.

Low-intensity light causes poor adhesion. Curing lights should be checked with reliable light testers at regular intervals. The emitting tip should be placed as close as possible to the curing surface.

- Place restorative material according to manufacturer's instructions.
- 2. Bonding of indirect restorations in combination with a light-curing luting cement: porcelain, Polyglas® and composite restorations (inlays, onlays, veneers, crowns).**
- 2.1. Pre-treatment of restoration**
- Pre-treat the bonding surface of the indirect restoration according to the manufacturer's instructions.

## 2.2. Dispensing and application of iBOND Self Etch

Light-cure the adhesive prior to application of the luting cement.

- Proceed as described in 1.
- Apply and process luting cement according the manufacturer's instructions. The luting cement must be completely light-cured after insertion of the indirect restoration.

## 3. Sealing hypersensitive areas of teeth.

### 3.1. Cleaning the tooth

- Clean tooth with an oil and fluoride-free polishing paste.

### 3.2. Dispensing and application of iBOND Self Etch

- Proceed as described in 1.
- Carefully remove oxygen-inhibited layer with an alcohol-soaked pellet after light-curing. If the desensitisation effect is not adequate, apply additional coats of iBOND Self Etch as described, light-cure and remove the oxygen-inhibited layer with a damp pellet.

### **Safety Instructions**

- The product is highly flammable.
- The product may irritate the eyes.
- In some cases an allergic reaction to ingredients of the product may occur in persons with known allergies.
- Contains methylacrylate and acetone.

### **Special Instructions**

- Keep away from children.
- The combination of iBond Self Etch with self-curing composites results in a significant reduction of adhesion and is therefore not recommended.
- Materials containing eugenol may inhibit the polymerisation of iBond Self Etch.
- Do not use iBond Self Etch after the expiry date.

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Heraeus Kulzer GmbH
Heraeus Kulzer, Inc.