DENTSPLY
DETREY

Prime&Bond[®] NT

Nano-Technology Dental Adhesive

Nano-Technologie Haftvermittler

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Prime&Bond[®] NT

Nano-Technology Dental Adhesive

Prime&Bond[®] NT is a universal self-priming dental adhesive designed to bond light-cured composite materials and Dyract[®] compomer materials to enamel and dentine as well as to metals and ceramic.

Prime&Bond NT combines primer and adhesive in a single bottle. The reduction of components and treatment steps simplifies use, maintaining superior bond strengths and protection against microleakage.

By mixing **Prime&Bond NT** with **Self-Cure Activator**, a universal, dual cure self-priming dental adhesive system is obtained. This system is designed to bond dual or self-cured composites to enamel, dentine and restorative materials even in areas inaccessible to light.

PRESENTATION FORM

Prime&Bond NT is available in a specially designed solvent-proof laminated bottle. Regular contents are 4.5 ml. The trial size bottle contains 1 ml only.

Alternatively, **Prime&Bond NT** is available pre-dosed for single patient use, each **quix[™]** - **unit dose** containing 0.125 ml.

COMPOSITION

- Di- and trimethacrylate resins
- Functionalised amorphous silica
- PENTA (dipentaerythritol penta acrylate monophosphate)
- Photoinitiators
- Stabilisers
- Cetylamine hydrofluoride
- Acetone

INDICATIONS

1. Self-priming bonding agent for:
 - composite restorations
 - Dyract compomer materials
 - cementation procedures for indirect restorations
 - adhesive repairs
2. Adhesive cavity varnish under amalgam restorations.
3. Protective varnish for hypersensitive cervical areas.

CONTRAINDICATIONS

Direct or indirect pulpcapping.

Use in patients with a known allergy to dimethacrylate resins.

WARNINGS

1. Prime&Bond NT contains methacrylates that may be irritating to the eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

2. Prime&Bond NT may cause sensitisation by contact with skin or mucous membranes in susceptible persons. After accidental contact, wash immediately with plenty of soap and water or rinse with plenty of water.
If sensitisation has occurred, discontinue use.
3. Prime&Bond NT contains acetone. Acetone is highly flammable. Keep away from sources of ignition – no smoking. Do not breathe vapour.
Take precautionary measures against static discharges.

PRECAUTIONS

Avoid Prime&Bond NT saturating gingival retraction cord. If Prime&Bond NT soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.

INTERACTIONS WITH DENTAL MATERIALS

Eugenol containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

If H₂O₂ has been used to clean the cavity, proper rinsing is essential. Higher concentration of H₂O₂ may interfere with the setting of polymerisable material and should not be used prior to the application of Prime&Bond NT.

Prolonged and intensive contact with acetone-containing products may lead to minute dissolution of the outermost surface of calcium hydroxide materials. This has no detrimental effect on the adhesion to the cavity walls.

ADVERSE REACTIONS

The following adverse reaction has been associated with the use of acetone solutions and acrylate monomers:

- Reversible inflammatory changes of the oral mucosa after accidental contact.

STEP-BY-STEP INSTRUCTIONS

1. Cleaning

Cavity cleanliness is paramount for the development of adhesion.

In cases where no cavity preparation has been made, clean the tooth surface with a rubber cup and pumice or a Prophy Paste like Nupro[®]. Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

2. Pulp Protection

For direct or indirect pulp-capping cover the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e.g. Dycal[®]), leaving the rest of the cavity surface free for bonding with Prime&Bond NT.

I. LIGHT-CURED COMPOSITE RESTORATIONS

3. Acid Conditioning of Enamel and Dentine (Total Etch Technique)¹

When used as a bonding agent for composite materials, it is recommended to follow the Total Etch Technique described below:

¹ Alternatively to the Total Etch Technique, the conventional Enamel Etch Technique can be followed. In this case, the enamel margins only are treated with DeTrey Conditioner 36 for at least 15 seconds. Then rinse and dry as described below.

3.1 Acid Conditioning

Apply 36% phosphoric acid (DeTrey® Conditioner 36) to the cavity surfaces starting at the enamel margins. For best results, condition enamel for at least 15 seconds and dentine for 15 seconds or less.

3.2 Rinsing and Drying

Remove gel with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.

Remove water from the rinsed cavity with a soft blow of air. Avoid desiccating the dentine, leave a moist surface.

Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with vigorous water-spray, dry, and repeat conditioning procedure of enamel for 5 seconds only. Rinse and dry as described above.

4. Application of Prime&Bond NT

1. Dispense Prime&Bond NT directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish or standard dappen dish.

For handling of quix - unit dose, please see illustrated instructions.

2. Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all cavity surfaces.

These surfaces should remain fully wet for 20 seconds and may necessitate additional application of mixed adhesive.

Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds.

Surface should have a uniform, glossy appearance.

If not, repeat application and air-dry.

3. Light-cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.

4. Immediately place composite restorative over the cured Prime&Bond NT.

II. DYRACT COMPOMER MATERIALS

3. Conditioning

For most restorative procedures with Dyract, it is not necessary to condition the prepared tooth. In this case, **proceed to step 4** of this section.

In the case of stress-bearing class I and II permanent as well as class IV Dyract®AP restorations and in situations where for cosmetic reasons bevelled enamel margins are preferred conditioning with NRC™ Non-Rinse Conditioner is recommended (see Directions for Use for NRC).

Alternatively to NRC, conditioning with phosphoric acid (DeTrey Conditioner 36) as described for composite restorations, may be applied.

4. Application of Prime&Bond NT

1. Dispense Prime&Bond NT directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish or standard dappen dish.

For handling of quix - unit dose, please see illustrated instructions.

2. Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all cavity surfaces.

These surfaces should remain fully wet for 20 seconds and may necessitate additional application of mixed adhesive.

Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds.

Surface should have a uniform, glossy appearance.

If not, repeat application and air-dry.

3. Light-cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.
4. Immediately place composite restorative over the cured Prime&Bond NT.

III. INDIRECT RESTORATIONS

1. Treatment of Restoration

Treat bonding surface of restoration according to manufacturer's or dental laboratory's instructions, e.g. etching or mechanical roughening followed by silanisation. Avoid touching the inner surface afterwards.

2. Cleaning and Conditioning of Preparation

Clean enamel and dentine prepared during a previous appointment with a rubber cup and pumice or a cleaning paste such as Nupro Prophylaxis Paste. Wash thoroughly with water spray and then air-dry. Clean freshly cut enamel and dentine with water spray and then air-dry. Condition dentine and enamel as described for direct restorations. In case of full crown preparations and the use of phosphoric acid, it is recommended to condition the enamel margins only.

3. Bonding with Prime&Bond NT

3.1 Application of Prime&Bond NT for use with a light-cured luting composite

3.1.1 Application to restoration

Apply Prime&Bond NT to the inner surface of the restoration with an Applicator Tip or disposable brush. Remove solvent by blowing with air from a dental syringe. Light-cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.

3.1.2 Application onto preparation

a) Dispense Prime&Bond NT directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish or standard dappen dish.

For handling of quix - unit dose, please see illustrated instructions.

b) Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all cavity surfaces.

These surfaces should remain fully wet for 20 seconds and may necessitate additional application of mixed adhesive.

Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance.

If not, repeat application and air-dry.

c) Light-cure for a minimum of 10 seconds². Ensure uniform exposure of all cavity surfaces.

d) Immediately place composite restorative over the cured Prime&Bond NT.

3.2 Application of Prime&Bond NT for use with a dual or self-cured luting composite

In these situations, Prime&Bond NT is mixed with Self-Cure Activator which is available as a separate product. For step-by-step instructions, please see directions for use for Self-Cure Activator.

3.3 Application of Prime&Bond NT for use with Dyract Cem^{plus}

For application of Prime&Bond NT follow the procedure described in the Directions for Use for Dyract Cem^{plus}.

² Check curing light for minimum curing output of at least 300 mW/cm².

IV. COMPOSITE REPAIRS

1. Roughen fractured composite surface as much as possible with a diamond bur and create mechanical retention in old composite, where possible.
2. For best results, microetch the surface of the old composite with an Intra-Oral Microetcher® (Danville Engineering) using 50 micron aluminium oxide powder to provide increased micromechanical retention. Rubber dam must be used throughout this procedure.
3. Wash thoroughly.
4. Air-dry with clean oil-free and moisture-free air.
5. Etch any exposed enamel margins and the surface of the composite with 36% phosphoric acid (DeTrey Conditioner 36) for 15 to 30 seconds. Rinse with water for 15 to 20 seconds, and air-dry.
6. Apply and light-cure Prime&Bond NT as described for direct restorations (see I.).
7. Complete repair with placement and curing of the appropriate shade of a suitable composite (e.g. Esthet•X®) or compomer (e.g. Dyract AP, Dyract® *flow*).

V. PORCELAIN/METAL REPAIRS

1. Some porcelain fractures on porcelain-fused-to-metal restorations do not expose metal, but most do.
2. Roughen the exposed metal and porcelain with a diamond instrument. Create mechanically retentive areas in the metal, if desired. For best results, sandblast the exposed metal and the porcelain margins with an Intra-Oral Microetcher (Danville Engineering) using 50 micron aluminium oxide to etch the surface of the metal. Rubber dam must be used throughout this procedure.
3. Bevel fractured porcelain margins with a clean diamond. Rinse and dry the preparations.
4. Apply silan (Calibra™ silane coupling agent) on the exposed porcelain margins following manufacturer's instructions.
5. Apply and light-cure Prime&Bond NT over the exposed porcelain and metal as described for direct restorations (see I.).
6. Complete repair with placement and curing of the appropriate shade of a suitable composite (e.g. Esthet•X) or compomer (e.g. Dyract AP, Dyract *flow*).

VI. AMALGAM REPAIRS

1. Roughen the old amalgam with a diamond instrument. Create mechanical retentive areas in the amalgam, if possible. For best results, sandblast the area to be repaired with an Intra-Oral Microetcher (Danville Engineering) using 50 micron aluminium oxide to etch the surface of the metal. Be sure to use rubber dam isolation.
2. Rinse microetched areas for 15 to 20 seconds with water, air-dry.
3. Etch any exposed enamel margins (but not the surface of the microetched amalgam) with 36% phosphoric acid (DeTrey Conditioner 36) for at least 15 seconds. Rinse with water for at least 15 seconds and air-dry.
4. Apply and light-cure Prime&Bond NT over the amalgam as described for direct restorations (see I.).
5. Complete repair with placement and curing of the appropriate shade of a suitable composite (e.g. Esthet•X) or compomer (e.g. Dyract AP, Dyract *flow*).

VII. ADHESIVE CAVITY VARNISH FOR USE WITH FRESH AMALGAM

1. Finish preparation.
2. If preparation is in close proximity to the pulp, place a hard-setting calcium hydroxide liner (e.g. Dycal).
3. Rinse and carefully air-dry cavity preparation, but do not desiccate exposed dentine.
4. Apply and light-cure Prime&Bond NT as described for direct restorations (see I.).
5. Place and condense amalgam (e.g. Dispersalloy®) as per usual technique.

Note: Acid etching of prepared cavity (enamel or enamel/dentine) is optional prior to placement of Prime&Bond NT when used as an adhesive cavity varnish.

VIII. COATING OF HYPERSENSITIVE CERVICAL AREAS

1. Clean cervical area of hypersensitive tooth with prophylaxis paste (Nupro) and rubber cup.
2. Remove prophylaxis paste with water spray. Dry cleaned area with a 2-second blow of air free of oil or water contamination. Avoid desiccating the dentine, leave a moist, but not wet glistening surface.
3. Apply and light-cure Prime&Bond NT as described for composite restorations (see I.).

STORAGE

Keep out of sunlight. To be stored at temperatures between 10 and 24 °C.
The Prime&Bond NT bottle should be tightly closed immediately after use.
Keep in a well ventilated place.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

Manufacturer:
DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE
Phone (0 19 32) 85 34 22

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© DENTSPLY DeTrey 2000-05-30