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# Exakto-Rock S

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## Instructions for use

# Exakto-Rock S

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

### 1. Description of the product

Exakto-Rock S is used for the fabrication of dental models with high precision and hardness. This super-hard stone can be processed with A-silicone, K-silicone, polyether, hydrocolloid and alginate impression materials.

This super-hard stone features high hardness, extended processing time and low setting expansion, which are specified in the requirements of DIN EN ISO 26873 for class IV stones. Due to its thixotropic modification, the excellent flowability provides it with a low viscosity on the vibrator.

High compressive strength and edge stability are preconditions for the required precision of fit. Special synthetic alpha plasters, which are produced in the continuous autoclaving process, are used as raw materials. Careful in-process controls ensure the high consistency of stone lots. The scannability of the stone enables processing in dental scanners and ensures perfect reproduction and recognition of all details.

### 2. Indication

Class IV stone for the fabrication of models, arches and master models. Also suitable for scanning.

### 3. Contraindication

Based on the current state of knowledge there are no known contraindications.

### 4. Safety recommendations and hazard warnings

For safe handling, please request the material safety sheet.

### 5. Storage and durability

#### 5.1 Storage

As far as storage of Exakto-Rock S is concerned, we recommend to ensure airtight sealing by double folding the bag.

#### 5.2 Durability

The durability of Exakto-Rock S is at least 2 years after date of manufacture if the bag is not opened.

### 6. Processing

- Mixing ratio: 20 ml distilled water / 100 g powder
- First fill water into the mixing bowl and then add plaster powder
- Let it soak at least 40 sec
- Premix manually for 20 sec
- Mix under vacuum for 40 to 60 sec
- Use vibrator for bubble-free and smooth flowing in
- Devest and process after 40 min

If alginates, hydrocolloids and gel duplicating materials that are suitable for plasters are used, the mixing time must be extended to 90 to 120 sec. Impressions do not have to be pretreated. After careful blow-drying, the surface of the impression should have a silky luster.

**bredent**

GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany · Tel. (+49) 0 73 09 / 8 72-22 Fax (+49) 0 73 09 / 8 72-24  
www.bredent.com · e-mail info@bredent.com

The quantity of water can be increased / reduced by 1.5 ml (for 100 g of powder). If less liquid is used, the stone will be harder and the processing time will be reduced. If more liquid is used, the stone will be „softer“ and the processing time is extended.

### 7. Technical data according to DIN EN ISO 26873

Mixing ratio	20 ml water / 100 g powder
Processing time	5 to 6 min at 18 to 20° C
Compressive strength	> 75 N/mm <sup>2</sup>
Setting expansion	< 0.08 %

### 8. Symbols

 Lot number

### 9. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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